

1161 Ringwood Ct. Suite#100, San Jose, California, 95131. U.S.A. Tel: 408-321-9898. 1-866-628-9898(USA) Fax: 408-321-9896.

E-mail: info@abbomax.com, www.abbomax.com

Ab^{Max} antibody affinity purification (Cat# 60005)

We provide the following purification according to customer requirements

- Peptide epitope affinity purification
- Protein A/G affinity purification
- MAP system affinity purification
- All the purified antibody will be verified by antigen (protein/peptide) specific ELISA
- Antibody enrichment and recovery efficiency

Sampsle	Prebleed		Test Bleed		AP Antibody		AP Antibody		Running Through		Protein A	
O.D	Mix of rabbit sera		Mix of Antisera		F1		F2		Waste		Purified Antibody	
Dilution	Test 1	Test 2	Test 1	Test 2	Test 1	Test 2	Test 1	Test 2	Test 1	Test 2	Test 1	Test 2
1:1K	0.309	0.312	1.596	1.609	1.908	1.919	1.08	1.075	0.659	0.667	1.114	1.124
1:5K	0.168	0.166	0.963	0.962	1.228	1.198	0.165	0.169	0.11	0.107	0.613	0.617
1:25k	0.061	0.06	0.274	0.277	0.687	0.685	0.069	0.066	0.084	0.077	0.177	0.167
1:125K	0.06	0.06	0.062	0.063	0.302	0.307	0.062	0.062	0.063	0.068	0.065	0.063
1:625K	0.06	0.058	0.059	0.059	0.106	0.104	0.061	0.064	0.074	0.065	0.066	0.06
1:3125K	0.057	0.06	0.059	0.06	0.061	0.063	0.061	0.06	0.063	0.066	0.063	0.062
1:15625K	0.058	0.06	0.061	0.058	0.062	0.067	0.064	0.068	0.067	0.062	0.061	0.064
1XPBS	0.056	0.059	0.059	0.058	0.064	0.065	0.064	0.064	0.062	0.065	0.069	0.063



